
Name:	Date	Section
Lab Partner: <u>ONLINE</u>	-	
<b>Experimental</b>	Observation	<u>s</u>
Part 1:		
T = <u>21.2</u> °C $P = $ <u>740.3</u>	torr V	= <u>0.500</u> L
T = K P =	atm	
Molar Mass of Butane ( $C_4H_{10}$ ) =		g/mole
Initial weight of cartridge; <u>67.</u>	<u>485</u>	g
final weight of cartridge: <u>66.</u>	352	_ g
mass of butane: g ;	n =	moles
Simplistic Calculations:		
1. apparent molar volume ( $V_m = V/n$ ) of butane as	t experimental T	$V_{\rm m} = \L/{\rm mole}$
2. apparent molar volume of butane at- STP:		$V_m = \_\ L/mole$
<b><u>Rigorous Calculations</u>:</b>		
3. Partial Pressure of water vapor in flask*:		$P_w = $ torr
4. Partial Pressure of butane in flask:		$P_{B} = P_{total} - P_{w} = $ torr
		= atm
5. Partial Pressure of butane (calc with vdW equa	tion**):	$P_{vdw} = $ atm
6. Compressibility Factor for Butane:		$Z_{\rm B} = P_{\rm B} V/n_{\rm B} RT = \_$
7. Estimated Second Virial Coefficient for Butane	e at room tempe	rature:
$Z_{B} = 1 + B_{B}n_{B}/V$ ;		$B_B = \_\ L/mole.$
*Equilibrium Vapor Pressure of Water: $\ln Pw(torr) = 20.94$	13 - 5300/T (20	- 30"C).

 $** P_{vdw} = nRT/(V - nb) - an^{2/}V^2 \quad ; \quad R = 0.08206 \quad L - atm/mol-K; \quad a = 14.47 \ atm-L^2/mole^2 \ , \quad b = 0.1226 \ L/mole = 0.122$ 

## Part 2. Buoyancy Effects:

(Note: This part of the experiment requires about 1 gram of butane.)

Record the initial weights of the deflated zip-lock bag (A, below) and the butane cartridge. In the fume hood, zip the bag as snugly as possible around the Tygon tubing (B). Connect the butane cartridge and inflate the bag to at least half, but no more than 80% of its capacity (C). Quickly remove the tubing and seal the bag (D). Record the final weights of the bag and the cartridge.



## Empty the Ziploc bag under the fume hood.

## Part 3. Conservation of Mass:

Weigh the ziploc bag (A, below). Determine if the weight is affected by whether or not the bag is inflated with air <u>(fill with air NOT your breath)</u>. Tare the bag, add **2.5 - 2.8 grams** of sodium bicarbonate (NaHCO<sub>3</sub>) to the bag (B), and record the mass. Seal **8 - 10 ml** of 50% acetic acid in snap-cap 10- ml vial. Put the vial in the bag, add about 10 ml of distilled water to dissolve the sodium bicarbonate (D), expel as much air *as* possible, seal, and record the mass of the bag



with its contents. Without unsealing the bag, open the vial and mix the acetic acid with the bicarbonate solution. Place the bag on the balance and observe the change in weight. After the effervescence has essentially ceased, record the weight of the bag and its contents.

Calculate the volume of gas expected from the reaction of sodium bicarbonate with excess acid.

mass of NaHCO<sub>3</sub> <u>2.685</u> g

molar mass of NaHCO<sub>3</sub> \_\_\_\_\_g/mole

moles of NaHCO<sub>3</sub> \_\_\_\_\_moles.

Weight of bag and reaction components:

before reaction:38.566gafter reaction:37.352gdiscrepancy:\_\_\_\_\_g

Using the density of air determined in Part 2, estimate the volume of gas produced by the reaction.

Estimated volume of expansion: \_\_\_\_\_ L (weight discrepancy/density of air)

## **Reaction:**

 $NaHCO_3(aq) + CH_3CO_2H(aq) -> _____ + ____ + _____ + _____$ 

Calculate the volume of gas expected from the reaction of sodium bicarbonate with excess acid.

Expected moles of CO<sub>2</sub>(g): \_\_\_\_\_ moles

Expected volume of gas at laboratory T & P: \_\_\_\_\_ L

 $\Delta V_{GAS} = (nCO_2RT/P_{ATM}) / [1 - (P_W/P_{ATM})]$